

Disposable Sack Sampler

Take Samples through the side of a Sack

The disposable Sack Sampler is a hygienic, single use device that has been specially designed to take samples through the side of a sack.

- Individually Bagged
- FDA & EU 10/2011 Conforming Materials
- BSE/TSE Free
- Full Batch Traceability
- ATEX Rated
- Available Sterile (gamma irradiated)

1. Insert the sampler through the side of the sack.
2. Angle sampler downwards to allow the product to flow into the sample collecting bag
3. Seal hole with an Ultra Adhesive label (see pages 138-143 for details)



QUOTE REQUEST ✓

BUY ONLINE ✓

| Stock No. | Length | Diameter | Material | Pre-Sterilized | Per Box |
|---------------|--|----------|----------|----------------|---------|
| 675 8085A100 | 400mm | 25mm | HDPE | No | 100 |
| 675 8085A100S | 400mm | 25mm | HDPE | Yes | 100 |
| 675 8085N100 | 400mm | 14mm | HDPE | No | 100 |
| 675 8085N100S | 400mm | 14mm | HDPE | Yes | 100 |
| 675 8080A450 | Sample Collecting Bag 450ml - for hygienic use - sterile | | | Qty: 1000 bags | |
| 675 8080A900 | Sample Collecting Bag 900ml - for hygienic use - sterile | | | Qty: 1000 bags | |