

Part No.


8050P-500

Description

Disposable LiquiThief – 500mm

Material of Construction:

Tip	PP (virgin, FDA approved – CFR 177.1520)
Outer Tube	PP (virgin, FDA approved – CFR 177.1520)
Inner Rod	PP (virgin, FDA approved – CFR 177.1520)
Clutch	PP (virgin, FDA approved – CFR 177.1520)
Piston Housing	PP (virgin, FDA approved – CFR 177.1520)
Piston Seal	Silicone Rubber (virgin, FDA approved – CFR 177.2600)
Handle	PP (virgin, FDA approved – CFR 177.1520)



Method of Construction:

Ultrasonic welding

Nominal Length

Fully closed: 523mm

Fully extended: 881mm

Diameter of Outer Tube:

21mm

Diameter of hole in tip:

1.5mm

Nominal Weight of Sampler:

79g

Assembly & Packing Environment:

Class 100,000 Medical Cleanroom

Individually Bagged?

Yes (heat sealed PE bag)

Method of Sterilisation:

N/A

Number of Samplers per Box:

20

Recommended Storage Conditions:

Dry and ambient temperature