


Disposable Powder Lance
Cross Sectional Sampler - No Cleaning Required

The disposable Powder Lance is ideal for taking a cross sectional sample through a bulk powder.

- Assembled & Packed in a Cleanroom
- Individually Bagged
- FDA & EU 10/2011 Conforming Materials
- BSE/TSE Free
- Full Batch Traceability
- ATEX Rated
- Available Sterile (gamma irradiated)



Without End Tip - ideal for use with sticky, cohesive powders



The disposable Powder Lance is quick and easy to use:

1. Push the disposable Powder Lance into the product, rotate the sampler and then withdraw
2. Deposit the sample into a suitable container
3. Place the used sampler into its plastic bag

Without End Tip

Stock No.	Pre-Sterile	Blade Length	Overall Length	Volume	Material	Per Box
675 8063H300	No	300mm	400mm	75ml	HDPE	20
675 8063H300S	Yes	300mm	400mm	75ml	HDPE	20
675 8063H500	No	500mm	600mm	125ml	HDPE	20
675 8063H500S	Yes	500mm	600mm	125ml	HDPE	20
675 8063H1000	No	850mm	950mm	200ml	HDPE	20
6758063H1000S	Yes	850mm	950mm	200ml	HDPE	20

QUOTE REQUEST 

BUY ONLINE 